BIOL 343 Applied Bioinformatics I

Post-alignment QC

Alignment is the most important step in RNA-seq analysis Counting (also important) and DEG ID relies on high-confidence mapping

- Recall the goal of our RNA-seq experiments...
 - Treatment vs Control
 - Mutant vs Wild type
 - Identify differentially expressed genes (DEGs)
- DEGs will be identified using statistical tests comparing *expression values* of transcripts/genes
- Expression values will be calculated based on the number of reads that *align/map* to a given genomic locus
- Filtering out spurious or redundant alignments is an important QC step after mapping


Browsing RNA-seq alignments JBrowse2

- 1. Open JBrowse2
- 2. Download the genome files (.fa & .gtf)
- 3. Download the STAR alignment (Aligned.sorted.bam) and indexes
- 4. Load the genome (.fa)
- 5. Add tracks:
 - 1. GTF
 - 2. Alignments

Browsing RNA-seq alignments JBrowse2

- 1. Navigate to SM_V10_4:12,110,496..12,137,215
- 2. Three dots on RNA-seq track...
 - 1. Pileup settings > Filter by...
 - 2. Pileup settings > Color by > Strand
 - 3. SNPCoverage settings > Draw arcs (deselect)

Read must have ALL th	nese flags	Read must have NONE of these flags
0		3588
□ read paired □ read mapped in proper pair □ read unmapped □ mate unmapped □ read reverse strand □ mate reverse strand □ first in pair □ second in pair □ not primary alignment □ read fails platform/vendor quality checks □ read is PCR or optical duplicate □ supplementary alignment		read paired read mapped in proper pair read unmapped mate unmapped read reverse strand mate reverse strand first in pair second in pair not primary alignment read fails platform/vendor quality checks read is PCR or optical duplicate supplementary alignment
	nple: filter tag na	e value field to get all reads containing any me SA with value * to get all ples include HP for haplotype, or RG for read
ŭ		



Browsing RNA-seq alignments JBrowse2

- 1. Close STAR Alignment
- 2. HISAT Alignment three dots...
 - 1. Pileup settings > Group by...
 - 1. Tag name = RG (read group)
 - 2. Deselect checkbox
- 3. Close the STAR alignment and all D/E alignments
- 4. Navigate to SM_V10_4:12,110,496..12,137,215
- 5. Three dots...
 - 1. SNPcoverage settings > remove arcs

